



Role of Textile Materials in Physical Medicine, Rehabilitation Studies and Reports

N Gokarneshan

Formerly Professor and Head, Department of Textile Chemistry, SSM College of Engineering, Komarapalayam, Tamil Nadu, India

Citation: N Gokarneshan (2026) Role of Textile Materials in Physical Medicine, Rehabilitation Studies and Reports. *J of Adv in Phy Med and Reh.* 1(3), 1-6. WMJ/JAPMR-105

***Corresponding author:** N Gokarneshan, Formerly Professor and Head, Department of Textile Chemistry, SSM College of Engineering, Komarapalayam, Tamil Nadu, India

Submitted: 24-07-2026

Accepted: 28-07-2026

Published: 28-08-2026

Textile materials are central to physical medicine and rehabilitation (PM&R), acting as dynamic interfaces for motion sensing, muscle stimulation, and therapeutic support. Advanced "smart" fabrics merge conductive polymers with traditional elastic materials to create non-invasive wearable sensors, allowing clinicians to track physical recovery and improve patient mobility in real-time.

Key Roles in Rehabilitation Studies and Reports

Gait & Motion Tracking: E-textiles utilize piezoresistive and conductive fibers (e.g., silver-coated nylon woven into spandex) to map joint movement and posture. These sensors provide vital kinetic and kinematic data, enabling studies to precisely monitor how patients recover motor functions after neurological events like strokes.

Compression Therapy: Specialized knitted textiles made of elastomeric yarns apply constant unit pressure on targeted body parts. Studies highlight their effectiveness in managing lymphedema, reducing swelling, and treating burn wounds by improving

sensory and motor

Electrostimulation & Pain Management: Conductive textile substrates are engineered into wearables capable of delivering Transcutaneous Electrical Nerve Stimulation (TENS). These are highly beneficial for patients experiencing hand spasticity or arthritis, actively stimulating muscles when traditional rigid electrodes cannot be used.

Thermotherapy: Bifunctional smart textiles combine motion monitoring with thermotherapy. By emitting controlled heat, they promote vasodilation and reduce muscle spasms.

Orthopedic Support: Knitted spacer fabrics are widely evaluated for stabilizing lighter joint injuries. They ensure proper moisture-wicking, ventilation, and pressure distribution, vastly improving skin condition and patient comfort during recovery.

Recent Findings and Evidence

Neurological Recovery: Research emphasizes how textile-based pressure sensors integrated into daily garments act as inexpensive, continuous monitoring terminals to analyze patient balance and fall risk.

Rheumatoid Arthritis: Recent clinical reports on biofunctional garments (embedded with materials like Germanium) and arthritis gloves (15–32 mmHg compression) show statistically significant reductions in joint pain, stiffness, and extracellular fluid retention.

Bi-Functional Textiles: Modern trials detail core-sheath sensing yarns (CSSYs) that are simultaneously able to harvest biomechanical energy and track recovery conditions.

Role of woven textiles

Woven textiles provide dimensional stability, strength, and structural support, making them foundational in physical medicine and rehabilitation (PM&R). Current studies highlight their critical role across compression therapy, tissue engineering, and sensory-integrated e-textiles, directly enhancing motor function recovery and remote patient monitoring.

Key Roles in Rehabilitation Studies and Reports

Compression Therapy & Orthopedics: Woven and elastomeric fabrics deliver programmed pressure, which is crucial for reducing swelling, managing varicose veins, and healing burn wounds. Studies show that these textiles significantly improve sensory-motor functions and reduce joint stiffness.

Tissue Engineering & Implants: Because of their superior dimensional stability, woven fabrics act as vital scaffolds for tissue regeneration and tendon or ligament repair. They are also used to reinforce hydrogels in regenerative medicine.

E-Textiles and Wearable Sensors: Recent reports detail the integration of conductive fibers into woven fabrics to monitor joint kinematics and gait retraining. By incorporating piezoresistive materials, these textiles can translate body movements into electrical resistance changes, providing therapists with real-time feedback on range of motion.

Self-Powered Therapeutics: Advanced woven structures utilizing triboelectric nanogenerators

(TENGs) and piezoelectric polymers are being studied to simultaneously harvest body movement energy and deliver targeted muscle stimulation.

Essential Characteristics and Advancements

To be utilized in physical medicine, woven textiles must meet strict clinical and material benchmarks:

Biocompatibility&Non-toxicity: Essential, especially for implantable woven structures (e.g., artificial tendons), to ensure the body does not reject the material.

Porosity & Moisture Management: Must allow for fluid absorption (wound care) or proper ventilation (orthopedic braces) to promote skin health.

Mechanical Strength: Must maintain their original form and tension even during rigorous joint movement or downstream chemical/thermal sterilization.

Role of knitted textiles

Knitted textiles are revolutionizing physical medicine by providing breathable, highly elastic, and seamlessly contouring "second-skin" platforms. In rehabilitation, they serve as active therapeutic garments, wearable sEMG (electromyography) electrode arrays, and dynamic sensors. They enhance recovery by tracking biomechanics and applying targeted neuromuscular stimulation.

Key Roles and Clinical Applications

Wearable Robotics & Actuation: Specialized knitted pneumatic fabrics integrate hollow channels that convert fiber contraction into controlled, dynamic clothing pressure (3–32 kPa). This acts as a soft-driven assist for weakened muscles.

Biomechanical Sensing: Digital knitting techniques like 3DKnITS use conductive yarns to create flexible strain and pressure sensors. These map joint kinematics, track gait, and evaluate force distribution during everyday mobility recovery.

Rehabilitation Garments: Customized knitted sleeves and tops use integrated textile electrodes to continuously measure upper limb muscle activity or deliver muscle stimulation. For example, studies highlight knitted gloves that safely apply traction resistance to fingers

to improve motor function and dexterity.

Compression & Scar Therapy: Because of their superior conformability, knitted spacer fabrics are critical for controlling edema and reducing the severity of post-burn scars by maintaining constant, programmatic pressure.

Patient Benefits

Seamless Comfort: Knitted AFO (Ankle-Foot Orthosis) interface socks, like SmartKnit, eliminate agonizing seams and bunching, which protects sensitive skin from friction and pressure ulcers.

Long-term Compliance: Transitioning from rigid, clinical equipment to smart, stretchable garments allows patients to conduct self-administered rehabilitation comfortably at home.

Hygiene and Breathability: The high porosity and open loops of knitted fabrics facilitate moisture-wicking and thermal regulation, ensuring the garments can be worn continuously without causing skin irritation.

Role of nonwoven textiles

Nonwoven textiles are vital in physical medicine and rehabilitation, offering hypoallergenic, breathable, and highly adaptable materials for patient care. They are critical for managing pressure ulcers, protecting surgical and orthopedic wounds, supporting hydrocolloid therapy, and providing comfortable, infection-resistant barrier support during physical therapies and mobility training.

Core Roles in Rehabilitation

Reports and studies highlight the following key functions of nonwoven materials in rehabilitation settings:

Advanced Wound and Scar Care: Nonwovens are the standard for non-implantable medical dressings. Their engineered porosity and moisture control capabilities actively promote healing. Innovations

like biodegradable hyaluronan nonwovens are reported to accelerate tissue regeneration and improve handling during topical applications.

Orthopedic and Therapeutic Padding: Nonwovens provide superior cushioning and friction reduction, protecting fragile skin, preventing pressure ulcers (bedsores), and supporting patients with skin conditions during prolonged immobility and physical medicine treatments.

Hygiene and Incontinence Management: Disposable nonwovens are heavily utilized in rehabilitation wards to maintain sanitary environments and protect sensitive skin from irritation caused by incontinence.

Surgical and Post-Surgical Support: Nonwoven fabrics are essential in personal protective equipment (PPE), surgical gowns, and post-surgery pads, ensuring 100% sterility and infection control. They also serve as internal scaffolds in regenerative medicine.

Key Material Benefits

Rehabilitation reports emphasize why these materials outperform traditional woven fabrics in clinical settings:

Engineered Consistency: Unlike woven fabrics which have weak spots where threads cross, nonwovens provide uniform barrier protection against pathogens and fluids.

Breathability & Comfort: Modern nonwovens often created using SMS (spunbond-meltblown-spunbond) techniques allow air and heat to circulate while blocking fluids and harmful microbes.

Hypoallergenic and Lint-Free: They minimize the risk of rashes, skin trauma, and wound contamination during physical manipulation or dressing changes.

Sterilization Adaptability: Nonwovens can easily withstand and maintain their structural integrity following sterilization via UV radiation, gas, or steam.